

Stat Strip[®] Xpress Operator Training

Glucose Monitoring System

Stat Strip[®]
Glucose Monitoring System



nova[®]
biomedical

Introduction

Welcome to the Nova Biomedical[®] StatStrip[®] Xpress Operator Training!

This will provide a comprehensive overview of the Stat Strip meter features and operating procedures including:

- System Overview
- Quality Control
- Patient Testing
- Additional procedures

Review questions will evaluate your comprehension and competency.

Training is typically completed in less than 30 minutes.

Let's Get Started!



Introduction

- Stat Strip Xpress is a portable handheld analyzer for whole blood glucose monitoring.
- Incorporates patented Multi-Well® system that corrects for hematocrit.
- Eliminates interferences from acetaminophen, ascorbic acid, maltose, galactose, xylose, and lactose.
- Provides accurate results regardless of sample's oxygen status.
- No calibration codes required.



Introduction


Product Features

- Results reported in 6 seconds
- Uses only 1.2 uL of sample
- END fill test strip
- Reportable range 10 – 600 mg/dL (0.6 - 33.3 mmol/L)
- On-Board Results Data Storage
- Large-font result display for easy viewing
- Long-life coin-cell battery



Introduction


The StatStrip Xpress Meter



Introduction

Battery

- Meter is powered by a single 3V coin-cell battery, 2450. The battery is housed at the rear of the meter behind the battery cover.
- Typical life of a new 3 V battery is about 600 patient &/or QC tests.
- Replace battery when indicated.



Glucose Test Strips

Always Date Vials When Opened

- 2 vials per box, 50 strips per vial.
- Keep vial tightly closed when not in use.
- Store strips in their original vial at room temperature.
- Expiration date is printed on the vial of test strips.
- Once opened the test strips are stable at room temperature for up to **180 days** or until the expiration date, whichever is first.
- Write the Date Opened and Discard Date on the vial.



Quality Control Solutions

Always Date Vials When Opened

- Three levels of Quality Control are available.
- Store controls at room temperature.
- Keep the vials tightly closed when not in use.
- Expiration date is printed on the control vials.
- Once opened, Quality Control is stable for up to **90 days** or until the expiration date, whichever is first.
- Write the Date Opened and Discard Date on the vial.



Using The StatStrip Meter

Quality Control Analysis




- To begin QC Analysis and to turn the meter ON, insert a Glucose Test Strip into the meter.
- Once Strip is inserted, the meter will show a Flashing Blood Droplet indicating it has turned on and ready for analysis.



Using The StatStrip Meter

Quality Control Analysis

- Press the Right or Left Arrow Buttons to scroll to the QC levels.
- Scroll to the correct QC level (number) as shown at bottom right of the screen.



Using The StatStrip[®] Meter

Quality Control Analysis

- Gently mix the control vial before each use and discard the first drop before sampling.
- Touch the end of the test strip to a drop of control solution until the test strip fills and the analysis countdown begins. An audible tone may be heard also.
- When applying Quality Control to the test strip, always make sure you hold the meter level or pointed slightly downward.
- Never stand the meter straight up when applying Quality Control material to the test strip.





6

Using The StatStrip Meter

Quality Control Analysis- Result

- The QC Result screen will display when analysis is complete (after 6 seconds).
- The result will automatically be stored in meter memory with a date/time identifier. Record result and verify it is within acceptable limits.
- Remove the test strip to turn the meter off.



Using The StatStrip Meter

Patient Analysis


- To begin Patient Analysis and to turn the meter ON, insert a Glucose Test Strip into the meter.
- Once Strip is inserted, the meter will show a Flashing Blood Droplet indicating it has turned on and ready for analysis.



Using The StatStrip Meter

Patient Analysis


- Clean the patient's finger according to hospital protocol.
- Dry thoroughly.
- Perform a finger stick according to hospital protocol.



Using The StatStrip Meter

Patient Analysis

- When the blood drop appears, touch the end of the test strip to the blood drop.
- The test strip is full when the countdown begins.
- When applying blood to the test strip, always make sure the test strip fills completely. Blood fills the strip from the end by capillary action.
- Make sure that you hold the meter level or pointed slightly downward. Never stand the meter straight up when applying blood to the test strip.
- Do not touch the strip to the blood drop a second time.
- Results appear in 6 seconds.





6

Using The StatStrip[®] Meter

Patient Analysis- Result



- The Patient Result screen will display when analysis is complete (after 6 seconds).
- The result will automatically be stored in meter memory with a date/time identifier. Record result manually.
- Remove the test strip to turn the meter off.



Using The StatStrip[®] Meter

Result- Above/Below Linear Range


- If result is low (LO- less than measurement range) or high (HI- greater than measurement range) the following message screens will appear.



Using The StatStrip Meter

Review Results

- The meter is capable of storing up to 400 total results (patient and QC). Results are stored and automatically deleted on a first-in, first-out basis.
- To view previously stored results and turn the meter on, press the MODE button once for less than 3 seconds.
- If there are no results in memory, the meter will display “- - -”.
- If there are results in memory, the most recent result will be displayed.
- Press the RIGHT or LEFT Arrow Buttons to scroll through results.
- The meter will display “END” upon reaching last stored result.
- Results are stored by date/time identifier.




Battery

Battery Indicator


Battery is Low:

- Flashing blood drop & battery icon shown when running samples.
- Battery charge sufficient for up to 10 tests. Change battery soon.



Battery is Very Low:


- Battery icon only shown at lower left of display when attempting to use meter.
- Meter will not operate properly and should not be used until battery is replaced.



Battery

Battery Replacement

- Turn the meter over.
- Remove the battery cover by lifting on the tab.
- Pull out the plastic tab below the battery cover which will also pull the battery out. Remove the battery.
- Replace with a new battery.
- NOTE: Be sure to place battery over the bottom section of the plastic tab/sleeve to assure easy removal the next time. Do place battery over top of plastic tab.
- Replace the battery cover. Meter software version will be display for 3 seconds after battery replacement.



Maintenance

Cleaning The Meter

- Clean the meter according to hospital protocol.
- Immediately follow with a water dampened cloth to remove all cleaning residue.
- Dry thoroughly with a soft cloth or lint-free tissue.
- Caution: Do NOT immerse the meter or hold the meter under running water.
- Caution: Do NOT spray the meter with a disinfectant solution.



Troubleshooting

Help Aids

Analysis Related Error Codes (displayed on meter):

- E3- Used Strip Error
 - The test strip was used previously or is defective. Repeat test with new test strip.
- E4- Short Sampler Error
 - Insufficient sample (control or blood) aspirated onto strip. Repeat test with new strip.
- E8- Bad Strip Error
 - The test strip is defective or bad. Repeat test with new test strip.
- E9- Bad Sample Error
 - Problem detected with sample. Repeat test with new test strip.



*This concludes the
StatStrip
Operator Training Program*